

Name: .....

Signature: ..... Index No. ....

545/3

**CHEMISTRY  
PRACTICAL**

**Paper 3**

July/August 2023

2 hours



**WESTERN JOINT MOCK EXAMINATIONS**

Uganda Certificate of Education

**CHEMISTRY PRACTICAL**

**Paper 3**

2 Hours

**INSTRUCTIONS TO CANDIDATES:**

Answer **both** questions.

Answers are to be written in the spaces provided in this booklet. Use **blue** or **black** ink ball pen only. Any work done in **pencil** will **not** be marked except drawings.

You are **not** allowed to use reference books. (i.e. text books, booklets on qualitative analysis etc).

**All working must be clearly shown.**

Mathematical tables and silent non-programmable calculators may be used.

For Examiners' Use Only		
	Marks	Examiners Initials
<b>Q.1</b>		
<b>Q.2</b>		
<b>Total</b>		

1. You are provided with the following:

**BA<sub>1</sub>**, which is a solution containing 12.3 grams per litre of the metal sulphate,  $YSO_4 \cdot nH_2O$ .


**BA<sub>2</sub>**, which is a solution of sodium hydroxide solution.

**BA<sub>3</sub>**, which is a 0.1M hydrochloric acid solution.

The metal sulphate reacts with sodium hydroxide solution according to the following equation:



Then insoluble hydroxide reacts with hydrochloric acid according to the following equation:



You are required to determine the number of moles of water of crystallization in the metal sulphate,  $YSO_4 \cdot nH_2O$ .

### Procedure:

Pipette 25.0 cm<sup>3</sup> (or 20.0 cm<sup>3</sup>) of **BA<sub>1</sub>** into a clean conical flask. Using a measuring cylinder, transfer an equal volume of **BA<sub>2</sub>** as in the procedure above into the conical flask containing **BA<sub>1</sub>**, shake the contents of the conical flask thoroughly for about 10 seconds and allow it to stand.

Label the resultant mixture **BA<sub>4</sub>**. Titrate the mixture **BA<sub>4</sub>** using solution **BA<sub>3</sub>** from the burette while shaking gently until the precipitate just dissolves to form a colourless solution.

Repeat the procedure until you obtain consistent results.

Record your results in the table below.

### RESULTS

Volume of pipette used..... cm<sup>3</sup> (½ mark)

Final burette reading (cm <sup>3</sup> )			
Initial burette reading (cm <sup>3</sup> )			
Volume of <b>BA<sub>3</sub></b> used (cm <sup>3</sup> )			

- (a) (i) State the volumes of **BA<sub>3</sub>** used to calculate the average volume. (7½ marks)  
(01 mark)

.....

- (ii) Calculate the average volume of **BA<sub>3</sub>** used. (2½ marks)

.....

.....

..... cm<sup>3</sup>

- (b) Calculate the number of moles of;

- (i) hydrochloric acid in **BA<sub>3</sub>** that reacted with the metal hydroxide in **BA<sub>4</sub>**.

(02 marks)

(ii) metal hydroxide in **BA<sub>4</sub>** that reacted with the hydrochloric acid in **BA<sub>3</sub>**. (02 marks)

(iii) metal sulphate that reacted with the sodium hydroxide in **BA<sub>2</sub>**. (02 marks)

(iv) metal sulphate (**YSO<sub>4</sub> · nH<sub>2</sub>O**) in one litre of solution. (02 marks)

(c) Determine the value of n in the formula **YSO<sub>4</sub> · nH<sub>2</sub>O**. (04 marks)  
(H = 1, O = 16, formula mass of YSO<sub>4</sub> = 120)

2. You are provided with substance **Y**, which contains **two** cations and **one** anion. Carry out the following tests to identify the cations and anion in **Y**. Identify any gas(es) that may be evolved. Record your observations and deductions in the table. (29 marks)

Tests	Observation	Deduction
(a) Heat <b>two</b> spatula end-fuls of <b>Y</b> strongly in a dry boiling tube and then allow it to cool.		
(b) To the residue obtained in (a), add 5 cm <sup>3</sup> of dilute nitric acid and shake well or warm to dissolve. Add dilute sodium hydroxide solution drop-wise until in excess. Filtrate and keep both the filtrate and the residue.		

<p>(c) Wash the residue from (b) and then add dilute nitric acid drop-wise until the residue just dissolves. Divide the solution into four parts and test as follows:</p> <p>(i) To the <b>first</b> part of the solution, add dilute sodium hydroxide solution drop-wise until in excess.</p>		
<p>(ii) To the <b>second</b> part of solution, add 4-5 drops of dilute sodium hydrogen carbonate solution and heat</p>		
<p>(iii) To the <b>third</b> part of solution, add 4-5 drops of dilute sodium sulphate solution.</p>		
<p>(iv) Use the <b>fourth</b> part to carry out a test of your own choice to confirm the cation in the residue. <b>TEST</b></p>		
<p>(d) To the filtrate obtained in (b), add dilute nitric acid drop-wise until the solution is just acidic. Divide the acidic solution into five parts:</p> <p>(i) To the <b>first</b> part of the acidic solution, add dilute sodium hydroxide solution drop-wise until in excess.</p>		
<p>(ii) To the <b>second</b> part of the acidic solution, add dilute ammonium hydroxide solution drop-wise until in excess.</p>		
<p>(iii) To the <b>third</b> part of the acidic solution, add 3-4 drops of dilute lead (II) nitrate solution.</p>		
<p>(iv) To the <b>fourth</b> part of the acidic solution, add dilute hydrochloric acid followed by excess barium chloride solution.</p>		
<p>(v) Use the <b>fifth</b> part of the acidic solution to carry out a test of your own choice to confirm the anion in <b>Y</b>.</p>		

- (e) (i) The cations in **Y** are.....
- (ii) The anion in **Y** is.....

**END**