

545/3  
UCE Chemistry  
Paper 3  
August, 2023  
2 hours



## UNNASE MOCK EXAMINATIONS

Name.....

### UGANDA CERTIFICATE OF EDUCATION

Uganda Certificate of Education

#### CHEMISTRY

PAPER 3

2 HOURS

#### INSTRUCTIONS TO CANDIDATES

- This paper consist of two questions
- All questions are compulsory. Answers are to be written in the spaces provided **ONLY**.
- You are not allowed to use any reference books (i.e text books, booklets on qualitative analysis etc)
- All working must be clearly shown
- Answers must be written using **blue** or **black ink only**

| For Examiner's use only |  |  |  |
|-------------------------|--|--|--|
| Q.1                     |  |  |  |
| Q.2                     |  |  |  |
| TOTAL                   |  |  |  |

1. You are provided with the following:-

BA1, which is a solution made by dissolving 2g of sodium hydroxide in 500cm<sup>3</sup> of water

BA2, which is a 0.05M solution of a strong acid W

You are required to determine the basicity of W by finding the moles of sodium hydroxide that reacted with one mole of W.

**Procedure:-**

Pipette 20 or 25cm<sup>3</sup> of BA1 into a clean conical flask. Then add 2 – 3 drops of phenolphthalein indicator and titrate the solution with solution BA2 from the burette until the end point. Repeat the titration 2 – 3 times to obtain consistent results. Enter your results in the table below:-

Results

**Volume of pipette used** .....cm<sup>3</sup>

( $\frac{1}{2}$  mark)

|  | 1 | 2 | 3 |
|--|---|---|---|
| Final burette reading (cm <sup>3</sup> )   |   |   |   |
| Initial burette reading (cm <sup>3</sup> ) |   |   |   |
| Volume of BA2 used (cm <sup>3</sup> )      |   |   |   |

(7  $\frac{1}{2}$  marks)

Titrate values to calculate the average volume of BA2 used

( $\frac{1}{2}$  mark)

.....  
Average volume of BA2 used cm<sup>3</sup>

(2  $\frac{1}{2}$  mark)

**Questions:-**

(a) Calculate the:-

(i) Molarity for BA1 ( $\text{Na} = 23$ ,  $\text{O} = 16$ ,  $\text{H} = 1$ )

(5 marks)

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

(ii) Number of moles of Sodium hydroxide in BA1 that reacted

.....

.....

.....

.....

.....

(iii) Number of moles of BA2 that reacted

(3 marks)

.....

.....

.....

.....

.....

(b) Determine the number of moles of sodium hydroxide that reacted with one mole of W

.....

.....

.....

.....

.....

.....

2. You are provided with substance M which contains two cations and one anion. Carry out the following test on M to identify the cations and anion. Identify any gases evolved. Record your observations and deductions in the table below:-

| TESTS   | OBSERVATIONS | DEDUCTIONS |
|---|--------------|------------|
| (a) Heat two spatula endfuls of M strongly until there is no further change   |              |            |
| (b) Dissolve two spatula endfuls of M in about 5cm <sup>3</sup> of water. Divide the resultant solution into six portions |              |            |
| (i) To the first portion add Sodium hydroxide solution drop wise until in excess and warm.                                |              |            |
| (ii) to the second portion, add ammonia solution drop wise until excess   |              |            |

|  |  |  |
|--|--|--|
| (iii) To the third portion, add 2 – 3 drops of sodium sulphate solution and warm   |  |  |
| (iv) Use the fourth position to carry out a test of your own to confirm one of the cations in M<br><br>Test:.....<br><br>..... |  |  |
| (v) To the fifth portion, add an equal volume of dilute nitric acid followed by 3 – 4 drops of lead(II) nitrate solution       |  |  |
| (vi) Use the sixth portion to carry out a test of your own to confirm the anion in M<br><br>Test:.....<br><br>.....            |  |  |

(c) (i) Cations in M ..... and.....

(ii) Anion in M.....and.....

**END**





## UNNASE MOCK EXAMINATIONS

UGANDA CERTIFICATE OF EDUCATION

Uganda Certificate of Education

CHEMISTRY

Chemistry Confidential S.4 UNNASE 2023

**Each candidate should be provided with:**

- 1 pipette (20 or 25cm<sup>3</sup>)
- 1 burette (50cm<sup>3</sup>)
- 2 conical flasks
- 2 beakers(250cm<sup>3</sup>)
- 6 test tubes
- 1 boiling tube
- 100cm<sup>3</sup> of BA1
- 60cm<sup>3</sup> of BA2
- 2g of M
- Phenolphthalein indicator
- Easy access to reagents for identifying cations and anions.
- Easy access to heat source
- BA1 is a solution made by dissolving 4g of sodium hydroxide in 1000cm<sup>3</sup> of water.
- BA2 is 0.15M hydrochloric acid.

- M is aluminum ammonium sulphate